



Your New Choice for Gluten Testing is Here!

Are you looking for a Gluten ELISA kit which is

- Based on **R5 antibody**, recommended by Codex Alimentarius?
- Approved by **AOAC** Research Institute Performance Tested MethodsSM Program and validated for **ISO 17025** accredited laboratories?
- **One of the most sensitive** test in the market with a great dynamic range (3-200 ppm of gluten)?
- Suitable for your **automation** needs?
- Including non-toxic **extraction buffer**?

Then your choice for gluten testing is
**SENSISpec INgezim Gluten R5
96 Well ELISA Kit**



SENSISpec INgezim Gluten R5

The **SENSISpec Ingezim Gluten R5 Test Kit** uses the R5 monoclonal antibody for specific detection of prolamins from wheat, rye and barley. Research studies have shown that the R5 antibody specifically targets the component of the gluten fraction that is believed to be directly responsible for celiac disease.

With its **limit of quantification of 3.4 ppm gluten (1.7 ppm gliadin)** the SENSISpec Ingezim Gluten R5 Test Kit is one of the most sensitive test kits in the market, creating additional safety for consumers with celiac disease.

The results of the validation study confirm the adequacy of the kit to be used as a reliable tool within an allergen risk management program and a **validated assay for ISO 17025 accredited laboratories**.

We are committed to support you with the implementation of the SENSISpec INgezim Gluten R5 ELISA kit in your laboratory.

For a trial offer and for more information please visit [this page](#).
Or contact your local Eurofins Technologies sales representative or distributor.